



Recombinant Leifsonia xyli subsp. xyli tRNA dimethylallyltransferase (miaA)

Product Code	CSB-EP721379LAAC
Abbreviation	miaA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AE04
Product Type	Recombinant Protein
Immunogen Species	Leifsonia xyli subsp. xyli (strain CTCB07)
Purity	≥85% (SDS-PAGE)
Sequence	MQLYRGMDIG TAKLPPGEWR GLPHHLFDVL DVTDEAAVAR YQPEARRVVQ EIRERGSTPI LVGGSGLYVS SVVDFRFPG TDTALRARLE AELAAQPGGT LFQRLARDP EAAKRIGSSN GRRIVRALEV AELTGAAVSG ALPEEPKPWW PVRVLGLAAP REELVQRLDA RVERMWAEGE LAEVEGLIPL GIERSVTARR AIGYAQALAE LTGELTRSQA QAQTQRLTRR YARRQLSWFK RYPGIHWLDY DDPQLVAAAL ARL
Source	E.coli
Target Names	miaA
Protein Names	Recommended name: tRNA dimethylallyltransferase EC= 2.5.1.75 Alternative name(s): Dimethylallyl diphosphate:tRNA dimethylallyltransferase Short name= DMAPP:tRNA dimethylallyltransferase Short name= DMATase Isopentenyl-diph
Expression Region	1-263
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.