



Recombinant *Desulfotalea psychrophila* Ribosome-binding factor A (rbfA)

Product Code	CSB-MP721412DAAA
Abbreviation	rbfA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AJY5
Product Type	Recombinant Protein
Immunogen Species	<i>Desulfotalea psychrophila</i> (strain LSv54 / DSM 12343)
Purity	≥85% (SDS-PAGE)
Sequence	MSVWDPKEIL DSIGLGKRTE VKRSTRVADA IRKELAVLLL EKIRDPKLQD VGISRVTATD DLTIVKIYYT IYGPHKHVKA VDKAFLAATG FIRSHLARTL NMRYTPNLQF CYDKTAEKVA EIENIFQEIA KERQDNGENS
Source	Mammalian cell
Target Names	rbfA
Protein Names	Recommended name: Ribosome-binding factor A
Expression Region	1-140
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.