



# Recombinant *Oryza sativa* subsp. *japonica* Sulfhydryl oxidase 1 (QSOX1)

<b>Product Code</b>	CSB-EP721446OFG-B
<b>Abbreviation</b>	QSOX1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6AUC6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	RSLGGREGPG EVDADAAVDL NATNFDAFLK ASLEPWAVVE FFAHWCPACR NYKPHYEKVA KLFNGRDAAH PGLILMARVD CASKVNIDLC NRFSVDHYPF LLWGPPTKFA SAKWDPKQEN NEIKLIDDGR TAERLLKWIN NQMKSSFSLE DKKYENENML PKNASDPEQI VQAIYDVEEA TAQALQIILE RKTIKPKNRD SLIRFLQILV ARHPSKRCRR GSAELLINF D DHWSSNLSLS SQEGSKLLES VAEENHWICG KEVPRGYWLF CRGSKSETRG FSCGLWVLMH SLTVRIGDGE SQSTFTSICD FIHNFFICEE CRKHFYEMCS SVSAPFRTAR ELSLWLWSTH NKVNMRLMKE EKDMGTGDPL FPKVTWPPNQ LCPSCYRSSK VTDGAVDWNE DAVYQFLVNY YGKKLVSSYK ETYMESLQQQ EKKIVSEDSS ISNAASVPIG AALGVAIASC TFGALACFWR AQQKNRKQRK NWN
<b>Source</b>	E.coli
<b>Target Names</b>	QSOX1
<b>Protein Names</b>	Recommended name: Sulfhydryl oxidase 1 EC= 1.8.3.2 Alternative name(s): Quiescin-sulfhydryl oxidase 1 Short name= OsQSOX1
<b>Expression Region</b>	31-513
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.