



Recombinant Rat Serine protease 40 (Prss40)

Product Code	CSB-MP721463RA
Abbreviation	Prss40
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AXZ6
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	QDYKPSQ TPPPTSDALL KPMGRVQKEI CGKTKFQGKI YGGSIKAER WPWQASLIFR GRHICGAVLI DKNWVASAAH CFKRSCLKPSD YRILLGYNEL SNPSNYSRQM TLSKVIVHED YNKLHSQEKD IVLIQLHLPV RYSSHIFPVC VPDQTTKEPS DESCWISGWG MVTDDKFLQA PFPLLDSEVF LMNDQECEAF FQTPQISITE YDAIKDDMIC AGDITNQKST CRGDSGGPLV CLLDSYWYLV GLASWSGACL EPIHSPSIFT RVSHFSDWIE KKKADTPDVD PSLAPLEETA PSLIGWRNYS AGTTLEPRIC TTLSSQVLL LQSIWLRIL
Source	Mammalian cell
Target Names	Prss40
Protein Names	Recommended name: Serine protease 40 EC= 3.4.21.- Alternative name(s): Testicular serine protease 2
Expression Region	34-369
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.