



# Recombinant Rat SH2 domain-containing protein 4A (Sh2d4a)

<b>Product Code</b>	CSB-YP721473RA
<b>Abbreviation</b>	Sh2d4a
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6AYC8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MLRQILSEMF IDPDLLAELS EEQKQILFYK MREEQIRRWK EREAAMERKE SLPVKSRPKK ENGKSVHWKL GADKQVWVWV MGEHHLDKPY DVLCDEILAE REHLRAAKDA ELRKAQSLQF TNSLKVKSLN CDLQAVKKTE PQNVTRKAAA EEATDQGPRAP ITPKKDDKAQ TKDLTKKKDS EELKQTEDEK TKQIYKNWKE DSEWQASLRK SKAADEKRRS LAKQAREDYK RLSQRGRSGD GLQNPLTGPQ KPKRPPLPPK PQFLQPLGSP PKSLGNQGGVI RTKTSSTQED IIRWFKEEQL PFRAGYQKNS DTIAPWFHGI LTLKKANELL STGMPGSFLI RVSEKIKGYA LSYLSEEGCK HFLIDASANS YSFLGVDQLQ HATLADLVEY HKEEPITSLG KELLFPCGQ QDKPPDYLEL FQ
<b>Source</b>	Yeast
<b>Target Names</b>	Sh2d4a
<b>Protein Names</b>	Recommended name: SH2 domain-containing protein 4A
<b>Expression Region</b>	1-422
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.