



Recombinant Gracilaria tenuistipitata var. liui Chromophore lyase CpcS/CpeS homolog (cpcS)

Product Code	CSB-MP721498GBAK
Abbreviation	cpcS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6B8Q1
Product Type	Recombinant Protein
Immunogen Species	Gracilaria tenuistipitata var. liui (Red alga)
Purity	≥85% (SDS-PAGE)
Sequence	MTYKNILDKL V GKWIHQRTS YFIHNQQIDY HQEEIKLKQI DNIYISTKDN NSLYQYELSN KINNQKIYYI FSKQEESEFG KLHKITNDEI KYRFTIHTH NCIKIESVRE NISYNEYIYL INNRFKITIS ILKRGQKYLA TSFISEIKIF N
Source	Mammalian cell
Target Names	cpcS
Protein Names	Recommended name: Chromophore lyase CpcS/CpeS homolog EC= 4.-.-.-
Expression Region	1-151
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.