



Recombinant *Debaryomyces hansenii* Ubiquitin-like-conjugating enzyme ATG10 (ATG10)

Product Code	CSB-YP721529DIS
Abbreviation	ATG10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6BJP0
Product Type	Recombinant Protein
Immunogen Species	<i>Debaryomyces hansenii</i> (strain ATCC 36239 / CBS 767 / JCM 1990 / NBRC 0083 / IGC 2968) (Yeast) (<i>Torulasporea hansenii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MIKDGKISNY QFNKHIDEFN ELLLRNLRLI PSYSKSGGSD ANLDIINNYK KELNSNTKNS KKTYILTRIC LSGSYLPNCL LIDFTLSYDD FYMVPVLYFR AFKDNSKSTS GNIDETRVTP IVSTEELVSN YYSVLGLSSD SNLGPTVTLD SHHLITDSSV WFYVHPCETL HRLREFMEAD NCLLPDSEQ LQVVKYLSIW YATYSLGGIF PSISLRPTSQ P
Source	Yeast
Target Names	ATG10
Protein Names	Recommended name: Ubiquitin-like-conjugating enzyme ATG10 EC= 6.3.2.- Alternative name(s): Autophagy-related protein 10
Expression Region	1-221
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.