



Recombinant *Debaryomyces hansenii* Pre-mRNA-splicing factor PRP46 (PRP46)

Product Code	CSB-BP721570DIS
Abbreviation	PRP46
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6BU94
Product Type	Recombinant Protein
Immunogen Species	<i>Debaryomyces hansenii</i> (strain ATCC 36239 / CBS 767 / JCM 1990 / NBRC 0083 / IGC 2968) (Yeast) (<i>Torulasporea hansenii</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MNFNIIKII CENIQLAPII NMVSYKELFN INDTDLIVPS DTLSEAVYDN SKLDYIFQNV NRRNRDEVED ESINIQTQLP KQQNQKISKD MSEVGSKQVS EIRPDSDSRW KLLRVMAGAH QGWVRSCTVD PVTNKWFVTG SSDSTIKIWD LASSNLKATI TGHIMGVRS L AVSSRYPYLF SGSEDKTVKC WDLERTNSSS GCQIRNYHGH VGGIYAMALH PELDLLFTGG RDSVIRVWDL RSRTEIMVLS GHRSDITSIA SQIGDPQIIT SSMDATIRLW DIRKATTQLA LTHHSKSIRS MAMHPQEMTM CSGDTSGLNK EWLLPGGELL NEFGHSGENK IINTLSINPS NNTLFSGYDD GRMEFYDYVS GDLLQSDATT PVTGSTESAI YASTFDMSG RLITCEGDKS IKIWGEE
Source	Baculovirus
Target Names	PRP46
Protein Names	Recommended name: Pre-mRNA-splicing factor PRP46 Alternative name(s): Pre-mRNA-processing protein 46
Expression Region	1-417
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.