



Recombinant *Debaryomyces hansenii* DNA replication complex GINS protein SLD5 (SLD5)

| | |
|--------------------------|---|
| Product Code | CSB-BP721586DIS |
| Abbreviation | SLD5 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q6BXX6 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Debaryomyces hansenii</i> (strain ATCC 36239 / CBS 767 / JCM 1990 / NBRC 0083 / IGC 2968) (Yeast) (<i>Torulasporea hansenii</i>) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MNKDLIDDDI LQEFDDSNVH SKSSKYGIRN SENIYDQLVA AMLNERMSPD VLPYKHELMK EVLTQLSNQQ QYLLDSHEYG DSNVESGIVT GDFKLQLMII ETDIERLNYL VRLYLRLRLA KIDKFTIHYI NETSNDDPTN DRLLSPEET EYMHKHFIL TQLYNNSFLK KMPHFLTLLD DTSGGQSMIS VPDINQPVFI KVITKVPIII NLDEDEDLEL VENGIYVVKY SLIKKYIEIG DIVLI |
| Source | Baculovirus |
| Target Names | SLD5 |
| Protein Names | Recommended name: DNA replication complex GINS protein SLD5 |
| Expression Region | 1-245 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |