



Recombinant *Yarrowia lipolytica* Pre-mRNA-splicing factor CWC2 (CWC2)

Product Code	CSB-MP721594YAP
Abbreviation	CWC2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6C007
Product Type	Recombinant Protein
Immunogen Species	<i>Yarrowia lipolytica</i> (strain CLIB 122 / E 150) (Yeast) (<i>Candida lipolytica</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MSEAPESIET SPVVTDPVAE SVQNQGVEVV QTEAPTEAQV VKTKSKKKRT KDPALYHTRK ARIQVDPDSV NTDDRPPQTG TVYNIWFNKW SGGDKEDEKF NQKKADGRCS IARDSGYTRA DKVPGSYFCL YFARGLCTQG HKCEFLHRLP VLTDMFSPTT DCFGRDRFFD YRDDMGGIGS ISRVNRTLYV GRIHVSDAAK AGALDEIVSR HFSEWGDVDR IRVLHDKGVA FVTYATEVNA QFAKEAMAHQ SLDSGEVLNV RWATQDPDPL AQAREQRRLE ENAAEAIKRL LPQEYVDELE GRAKSKPLL EGYEEDDSEK LKRIMNNQKE AGAEPAAEEVK QIEPAPEPAA PVSTDTGMFN KSSLSALKAL KKKKKTAKV EE
Source	Mammalian cell
Target Names	CWC2
Protein Names	Recommended name: Pre-mRNA-splicing factor CWC2
Expression Region	1-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.