



Recombinant *Yarrowia lipolytica* Nascent polypeptide-associated complex subunit beta (EGD1)

Product Code	CSB-MP721603YAP
Abbreviation	EGD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6C2F3
Product Type	Recombinant Protein
Immunogen Species	<i>Yarrowia lipolytica</i> (strain CLIB 122 / E 150) (Yeast) (<i>Candida lipolytica</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MPVDPAKLAA LQKKSGAQSG GKGTPRRPGK KVAGRNISED EKKLSATLKK FNAQEITGIS EVNMFKEDGT VLHFPKVHVE GSVASNTFAI SGPSQQKDIA ELIPDILPQM GQDALLQLQQ AAVQFSKLQE QAKKTAGGPD AAKEAGDDEI PNLVENFEDN VE
Source	Mammalian cell
Target Names	EGD1
Protein Names	Recommended name: Nascent polypeptide-associated complex subunit beta Short name= NAC-beta Alternative name(s): Beta-NAC
Expression Region	1-162
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.