



Recombinant Kluyveromyces lactis Mitochondrial morphogenesis protein SLD7 (SLD7)

Product Code	CSB-BP721667KBK
Abbreviation	SLD7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6CJS0
Product Type	Recombinant Protein
Immunogen Species	Kluyveromyces lactis (strain ATCC 8585 / CBS 2359 / DSM 70799 / NBRC 1267 / NRRL Y-1140 / WM37) (Yeast) (Candida sphaerica)
Purity	≥85% (SDS-PAGE)
Sequence	MTLEKELILR LNVGHNTVIR DVQLWREFED GDGQIAVKNC LAEVVGMVNI DKLPFWLNRG QRYQCFTTSE TSVNYFRAKL MRQRHRNRGI VCKVTGSQVL TKEQFEYFLF YKRLKTDINA SFEIESMVLH LSTKHELDKK LIPLSADELQ SVVPSNAGAT MSIESILEKN KRQRILKSNS LNKQVLLNDQ QQQFMSLLSQ CILSGLRLRG VHRTQHEKLY KMTYKASEFA FRDELRETEP ISFESIQDCV ETLLKLFTKT
Source	Baculovirus
Target Names	SLD7
Protein Names	Recommended name: Mitochondrial morphogenesis protein SLD7
Expression Region	1-260
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.