



# Recombinant *Kluyveromyces lactis* Ribonuclease T2-like (RNY1)

<b>Product Code</b>	CSB-EP721707KBK-B
<b>Abbreviation</b>	RNY1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6CRT6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Kluyveromyces lactis</i> (strain ATCC 8585 / CBS 2359 / DSM 70799 / NBRC 1267 / NRRL Y-1140 / WM37) (Yeast) ( <i>Candida sphaerica</i> )
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	FD LQHIFYKEQY SCPISLPASC SNNTEIKDSC CFEFPGGIML QTQFWDYIPP KGVSDPDELV RHLGPLDSFT NHGLWPDNCD GTYAQFCNRE SNIDDVHLL NDDQFNRRDD LPINGTDLLE TMDMYWKSNT GDDESLWVHE YNKHGTICRT LYPCYKKGW VAGNSKKQAV YDYFRIAMKL FHKDQTYQTL KSAGIEPSVE KSYTKLEISN ALKEGHSGEV VHFLCDRHGS LNQIWFHSL KGSLLGEKFV PISALPKGSN CPDDNIKWYP KGHVPSSYRP PNGNHPGTRG VVRIGNGKGF LIKNGHWYLK GTPANFFLIE APFGNYLKT RMGYCGIDNS NKVLACNKNV AQAAQFEYDA KKGYIGYNGA YDWYATKYPR GNQQAPVYAG SNDDGYNFQL KFKVA
<b>Source</b>	E.coli
<b>Target Names</b>	RNY1
<b>Protein Names</b>	Recommended name: Ribonuclease T2-like Short name= RNase T2-like EC= 3.1.27.1
<b>Expression Region</b>	19-425
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.