



Recombinant Kluyveromyces lactis DASH complex subunit DUO1 (DUO1)

Product Code	CSB-YP721724KBK
Abbreviation	DUO1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6CUD4
Product Type	Recombinant Protein
Immunogen Species	Kluyveromyces lactis (strain ATCC 8585 / CBS 2359 / DSM 70799 / NBRC 1267 / NRRL Y-1140 / WM37) (Yeast) (Candida sphaerica)
Purity	≥85% (SDS-PAGE)
Sequence	MVEHKLDAST ISKMIPQIFD QMRRQTGSKN FASTPASIST QSLLKEQEQL DKIIPVIQNL NERIKDCKEG DMERIKETCH AMNTILDKFI SIQSQASYVN SMMNNEEYLD YVASGKSEEE YISAKQREVE ELEKRVKQYT TMHAPSVKPP ASKIPVGNDT RKVTKPGQRN GRSGLYSIST RNNSNNSGAN GRRPPGNGNN VGRRRVYR
Source	Yeast
Target Names	DUO1
Protein Names	Recommended name: DASH complex subunit DUO1 Alternative name(s): Outer kinetochore protein DUO1
Expression Region	1-208
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.