



Recombinant *Erwinia carotovora* subsp. atroseptica UPF0304 protein ECA3037 (ECA3037)

Product Code	CSB-YP721784EAAB
Abbreviation	ECA3037
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6D2Q8
Product Type	Recombinant Protein
Immunogen Species	Pectobacterium atrosepticum (strain SCRI 1043 / ATCC BAA-672) (<i>Erwinia carotovora</i> subsp. atroseptica)
Purity	≥85% (SDS-PAGE)
Sequence	MEMTNAQRLLSNQYKMMTM LDPDNAERYR RLQTIIEGY GLQMRELDRE FGELTEEVCR TIINVMEMHH ALQVSWTNLK DKQDLDERRL TFLGFDAATE ARYLGFVRFM VHVEGRYPHF DAGTHGFNAQ TKMWEKYTRM LAVWQSCPRQ YHLSAVEIAQ IINA
Source	Yeast
Target Names	ECA3037
Protein Names	Recommended name: UPF0304 protein ECA3037
Expression Region	1-164
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.