



Recombinant *Xenopus tropicalis* Nuclear cap-binding protein subunit 2 (ncbp2)

Product Code	CSB-BP721865XBF
Abbreviation	ncbp2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DES0
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus tropicalis</i> (Western clawed frog) (<i>Silurana tropicalis</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MALGALKSDS YVELSQYRDQ HFRGNRSDQE CLLKHSTLY VGNLSFYTTE EQIHELFSKS GDVKKIVMGL DKIKKTACGF CFVEYYTRTD AEQAMRFING TRLDDRIIRT DWDAGFKEGR QYGRGKSGGQ VRDEYRQDYD AGRGGYGKIV QSI
Source	Baculovirus
Target Names	ncbp2
Protein Names	Recommended name: Nuclear cap-binding protein subunit 2 Alternative name(s): 20 kDa nuclear cap-binding protein NCBP 20 kDa subunit Short name= CBP20
Expression Region	1-153
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.