



Recombinant *Xenopus tropicalis* Ragulator complex protein LAMTOR2 (lamtor2)

Product Code	CSB-BP721868XBF
Abbreviation	lamtor2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DF40
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus tropicalis</i> (Western clawed frog) (<i>Silurana tropicalis</i>)
Purity	>85% (SDS-PAGE)
Sequence	MLRPKALTQV LSQANTGGVQ STLLLNNEGS LLAYSGYGDT DARVTAAIAS NIWAAYDKNG HQAFNEDNLK FILMDCMEGR VAITRVSPLL LCMYAKETVG FGMLKAKAQA LVYYLEESLN QVSSS
Source	Baculovirus
Target Names	lamtor2
Protein Names	Recommended name: Ragulator complex protein LAMTOR2 Alternative name(s): Late endosomal/lysosomal adaptor and MAPK and MTOR activator 2 Roadblock domain-containing protein 3
Expression Region	1-125
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.