



# Recombinant Mouse Protein FAM5B (Fam5b)

<b>Product Code</b>	CSB-BP721875MO
<b>Abbreviation</b>	Fam5b
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6DFY8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	VSATVAA VVPEQHVSSA GQAPLDWLLT DRGPFHRAQE YADFMERYRQ GFTTRYRIYR EFARWKVNNL ALERRDFFSL PLPLAPEFVR NIRLLGRRPN LQQVTENLIK KYGTHFLLSA TLGGEESLTI FVDKRKLSRK SETTGGIPVV GGTGNSSAVS LETLHQLAAS YFIDRESTLR RLHHIQIATG AIKVTETRTG PLGCSNYDNL DSVSSVLVQS PENKVQLLGL QVLLPEHLRE RFVAAALSYI TCSSEGELVC RENDCWCKCS PTFPECNCPD ADIQAMEDSL LQIQDSWATH NRQFESEEF QTLLKRLPSD RFLNSTAISQ YWTMDSNLQH RYQQLGASLK VLLKKMHRIV RRLFNLCKRC HRQPRFRLPK ERSLSFWWNR IQSLLYCGES TFPGTFLEQS HSCTCPYDQS SCQGPIPCAL GEGPACAHCA SDNSTRCGSC NPGYVLAQGL CRPEVAESLE NFLGLETDLQ DLELKYLLQK RDSRIEVHSI FISNDMRLGS WFDPSWRKRM LLTLKSNKYK PGLVHVMLAL SLQICLTKNS TLEPVMAIYV NPFGGSHSES WFMPVNEGNF PDWERTNVDA AAQCQNWTIT LGNRWKTFE TVHVYLRSLR KSLDDSSNET IYYEPLMTD PSKNLGYMKI NTLQVFGYSL PFDPAIRD LILQLDYPYQ GSQDSALLQL IELRDRVNQL SPPGKVRLDL FSCLLRHRLK LANNEVGRIQ SSLRAFNSKL PNPVEYETGK LCS
<b>Source</b>	Baculovirus
<b>Target Names</b>	Brinp2
<b>Protein Names</b>	Recommended name: Protein FAM5B
<b>Expression Region</b>	34-783
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.