



Recombinant Mouse Protein FAM5B (Fam5b)

Product Code	CSB-MP721875MO
Abbreviation	Fam5b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DFY8
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	VSATVAA VVPEQHVSSA GQAPLDWLLT DRGPFHRAQE YADFMERYRQ GFTTRYRIYR EFARWKVNNL ALERRDFFSL PLPLAPEFVR NIRLLGRRPN LQQVTENLIK KYGTHFLLSA TLGGEESLTI FVDKRKLSRK SETTGGIPVV GGTGNSSAVS LETLHQLAAS YFIDRESTLR RLHHIQIATG AIKVTETRTG PLGCSNYDNL DSVSSVLVQS PENKVQLLGL QVLLPEHLRE RFVAAALS YI TCSSEGELVC RENDCWCKCS PTFPECNCPD ADIQAMEDSL LQIQDSWATH NRQFESEEF QTLLKRLPSD RFLNSTAISQ YWTMDSNLQH RYQQLGASLK VLLKKMHRIV RRLFNLCKRC HRQPRFRLPK ERSLSFWWNR IQSLLYCGES TFPGTFLEQS HSCTCPYDQS SCQGPIPCAL GEGPACAHCA SDNSTRCGSC NPGYVLAQGL CRPEVAESLE NFLGLETDLQ DLELKYLLQK RDSRIEVHSI FISNDMRLGS WFDPSWRKRM LLTLKSNKYK PGLVHVMLAL SLQICLTKNS TLEPVMAIYV NPFGGSHSES WFMPVNEGNF PDWERTNVDA AAQCQNWTIT LGNRWKTFE TVHVYLRSLR KSLDDSSNET IYYPELEMTD PSKNLGYMKI NTLQVFGYSL PFDPAIRDIL ILQLDYPYDQ GSQDSALLQL IELRDRVNQL SPPGKVRDL FSCLLRHRLK LANNEVGRIQ SSLRAFNSKL PNPVEYETGK LCS
Source	Mammalian cell
Target Names	Brinp2
Protein Names	Recommended name: Protein FAM5B
Expression Region	34-783
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.