



Recombinant Human Protein CC2D2B (CC2D2B)

Product Code	CSB-BP721888HU
Abbreviation	CC2D2B
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q6DHV5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MMTEKHEDHC LKSCSGHSYI RKNWLGCVF PFSALLQQSE FLDQTEVLQR AQIFKKNCKA MFPNRRIVTT VFNDEGIQFL VTRYIKALNP PQLLDIFLH NSNATFDLIA RFVSLIPFVP NTPDENDGSD IWMTSEHCIS LAIGNKEEHA ILLCNFFLYF GKKALVLLGT SVLEGHVAYV VTQETNEYLL WNPSTGQCYK QFDFPCPLKS VDCLFDDRNV WFNIQQNNTM MAVFFDYSKE SFWKQLLPKN VQGTKIQSIQ VTGFPIQMPY IDVQSIIDAV YQTGIHSAEF PQTEFALAVY IHPYPNNILS VVVYLASLVQ HQ
Source	Baculovirus
Target Names	CC2D2B
Protein Names	Recommended name: Protein CC2D2B
Expression Region	1-322
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>