



Recombinant Mouse F-box only protein 27 (Fbxo27)

Product Code	CSB-BP721890MO
Abbreviation	Fbxo27
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DIA9
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MGAWISRTRV PTPEPDPQEV LDLSRLPPEL LLLVLSHVPP RTLLMHCRRV CRAWRALVDG QALWLLLLAR DHSAAGRALL TLARRCLPPA HEDTPCPLGQ FCALRPLGRN LISNPCGQEG LRKWMVRHGG DGWVVEKNRK PVPGAPSQTC FVTSFSWCRK KQVVDLVEKG LWPELLDSGG VEIAVSDWWG ARHDSGCKYR LFVTLLEDAHQ NVIDKFSAVP DPIEQWNNDI YLQVTHVFSG IRRGIRFVSF EHWGQDTQFW AGHYGARVTN SSVIIRVCQS
Source	Baculovirus
Target Names	Fbxo27
Protein Names	Recommended name: F-box only protein 27
Expression Region	1-280
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.