



Recombinant *Mesoplasma florum* Ribosomal RNA small subunit methyltransferase G (rsmG)

Product Code	CSB-MP721986MEO
Abbreviation	rsmG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6F0F5
Product Type	Recombinant Protein
Immunogen Species	<i>Mesoplasma florum</i> (strain ATCC 33453 / NBRC 100688 / NCTC 11704 / L1) (<i>Acholeplasma florum</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MWTMFNNWKI FEDNLGFVPS EEIKDKLNEY YKILVEENLK YNLTRITNEE DVYEKHFLLDS LLFTKETKID NQKIIDIGTG PGFPGIVLKI FFPETDITLI DSNNKKINFL NIVIKKLNK QIEAKHARAE ELARIENEKY DIAISRAVAY LDVILELAVR FLKIQGKLIL LKGPRADEEI KNSKNIDQKL KIKLTNKQSL ADTGFGERIN LFYEKEKSTP ELYPREYAKI VKESGKK
Source	Mammalian cell
Target Names	rsmG
Protein Names	Recommended name: Ribosomal RNA small subunit methyltransferase G EC= 2.1.1.- Alternative name(s): 16S rRNA 7-methylguanosine methyltransferase Short name= 16S rRNA m7G methyltransferase
Expression Region	1-237
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.