



Recombinant Anaplasma marginale DNA-directed RNA polymerase subunit omega (rpoZ)

Product Code	CSB-YP722030AAAC
Abbreviation	rpoZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5PB40
Product Type	Recombinant Protein
Immunogen Species	Anaplasma marginale (strain St. Maries)
Purity	≥85% (SDS-PAGE)
Sequence	MREGDSGYFS CYDGNRFR LV ILASQRAHEL SSGACTAVAR KGDKNTVVAL REIVGEQLDL AAVFRLAVNR CRKYLEEFTN AREVAAARSS QAAPKSAPGQ EIGKSFREKD PSAAAF LDQE QFFSGGGGE
Source	Yeast
Target Names	rpoZ
Protein Names	Recommended name: DNA-directed RNA polymerase subunit omega Short name= RNAP omega subunit EC= 2.7.7.6 Alternative name(s): RNA polymerase omega subunit Transcriptase subunit omega
Expression Region	1-128
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.