



Recombinant Salmonella paratyphi A Ribosome maturation factor RimM (rimM)

Product Code	CSB-EP722069SBM
Abbreviation	rimM
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5PFF7
Product Type	Recombinant Protein
Immunogen Species	Salmonella paratyphi A (strain ATCC 9150 / SARB42)
Purity	≥85% (SDS-PAGE)
Sequence	MSKQLAAQVP AEPVVLGKMG SSYGIRGWLR VFSSTEDAES IFDYQPWFIQ KAGWQQVQL ESWKHNNQDL IIKLKGVDDR DAANLLTNCE IVVDSSQLPA LEEGDYWWD LMGCQVVTAE GYDLGKVIDM METGSNDVLV IKANLKDAFG IKERLVPFLD GQVIKKVDLA TRTIEVDWDP GF
Source	E.coli
Target Names	rimM
Protein Names	Recommended name: Ribosome maturation factor RimM
Expression Region	1-182
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.