



Recombinant Salmonella paratyphi A Ribosomal RNA large subunit methyltransferase I (rlmI)

Product Code	CSB-BP722076SBM
Abbreviation	rlmI
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5PGB6
Product Type	Recombinant Protein
Immunogen Species	Salmonella paratyphi A (strain ATCC 9150 / SARB42)
Purity	≥85% (SDS-PAGE)
Sequence	MTESTFPQYP RLVLSKGREK SLLRRHPWVF SGAVSRLEGK ANLGETIDIV DHQGKWLARG AWSPASQIRA RVWTFDKAES IDIAFFTRRL RQAQQWRDWL AKKDGLDSYR LIAGESDGLP GVTIDRFGYF LVLQLLSAGA EYQRAALISA LQTCYPDCAI YDRSDVAVRK KEGMALTQGP VTGELPPALL PIEEHGMKLL VDIQGGHKTG YYLDQRDSRL ATRRYVENQR VLNCFSTYGG FAVSALMGGC RQVVSVDTSQ DALDIARQNV ELNQLDLSKA EFVRDDVFKL LRAYREHGEK FDVIIMDPPK FVENKSQLMG ACRGYKDINM LAIQLLNPGG ILLTFSCSGL MTSDLFQKII ADAAIDAGR D VQFIEQFRQA ADHPVIATYP EGLYLKGFAC RVM
Source	Baculovirus
Target Names	rlmI
Protein Names	Recommended name: Ribosomal RNA large subunit methyltransferase I EC=2.1.1.191 Alternative name(s): 23S rRNA m5C1962 methyltransferase rRNA (cytosine-C(5)-)-methyltransferase RlmI
Expression Region	1-403
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.