



# Recombinant Salmonella paratyphi A Ribosomal RNA large subunit methyltransferase I (rlmI)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-MP722076SBM                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Abbreviation</b>      | rlmI                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                 |
| <b>Uniprot No.</b>       | Q5PGB6                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Immunogen Species</b> | Salmonella paratyphi A (strain ATCC 9150 / SARB42)                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Sequence</b>          | MTESTFPQYP RLVLSKGREK SLLRRHPWVF SGAVSRLEGK ANLGETIDIV<br>DHQGKWLARG AWSPASQIRA RVWTFDKAES IDIAFFTRRL RQAQQWRDWL<br>AKKDGLDSYR LIAGESDGLP GVTIDRFGYF LVLQLLSAGA EYQRAALISA<br>LQTCYPDCAI YDRSDVAVRK KEGMALTQGP VTGELPPALL PIEEHGMKLL<br>VDIQGGHKTG YYLDQRDSRL ATRRYVENQR VLNCFSTYGG FAVSALMGGC<br>RQVVSVDTSQ DALDIARQNV ELNQLDLSKA EFVRDDVFKL LRAYREHGEK<br>FDVIIMDPPK FVENKSQLMG ACRGYKDINM LAIQLLNPGG ILLTFSCSGL<br>MTSDLFQKII ADAIDAGRQ VQFIEQFRQA ADHPVIATYP EGLYLKGFAC<br>RVM |
| <b>Source</b>            | Mammalian cell                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Target Names</b>      | rlmI                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Names</b>     | Recommended name: Ribosomal RNA large subunit methyltransferase I EC=2.1.1.191 Alternative name(s): 23S rRNA m5C1962 methyltransferase rRNA (cytosine-C(5)-)-methyltransferase RlmI                                                                                                                                                                                                                                                                                                |
| <b>Expression Region</b> | 1-403                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                        |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life                                                                                                                                                                                                                                                  |



of lyophilized form is 12 months at -20°C/-80°C.