



Recombinant Salmonella paratyphi A Ribonuclease 3 (rnc)

Product Code	CSB-EP722095SBM-B
Abbreviation	rnc
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5PIG5
Product Type	Recombinant Protein
Immunogen Species	Salmonella paratyphi A (strain ATCC 9150 / SARB42)
Purity	≥85% (SDS-PAGE)
Sequence	MNPIVINRLQ RKLGYTFNHQ ELLQQALTHR SASSKHNERL EFLGDSILSF VIANALYHRF PRVDEGDMSR MRATLVRGNT LAELAREFDL GECLRLGPGE LKSGGFRRES ILADTVEALI GGVFLDSNIQ TVEQLILNWX KTRLDEISPG DKQKDPKTRL QEYLQGRHLP LPSYLVVQVR GEAHQDEFTI HCQVSGLSEP VVG TGSSRRK AEQAAA EQAL KKLELE
Source	E.coli
Target Names	rnc
Protein Names	Recommended name: Ribonuclease 3 EC= 3.1.26.3 Alternative name(s): Ribonuclease III Short name= RNase III
Expression Region	1-226
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.