



# Recombinant Salmonella paratyphi A Ribose import ATP-binding protein RbsA (rbsA), partial

<b>Product Code</b>	CSB-EP722115SBM
<b>Abbreviation</b>	rbsA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5PJX5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Salmonella paratyphi A (strain ATCC 9150 / SARB42)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDALLQLKGI DKAFPGVKAL SGAALNVYPG RVMALMGENG AGKSTMMKVL TGIYTRDAGS LLWLGKETTFF NGPKSSQEAG IGIIHQELNL IPQLTIAENI FLGREFVNRF GKIDWKKMYA EADQLLAKLN LRFKSDKLVG ELSIGDQQMV EIAKVLSEFES KVIIMDEPTD ALDTEETESL FRVIRELKSQ RRGIVYISHR MKEIFEICDD VTVFRDGQFI AEREVATLTE DSLIEMMVGR KLEDQYPHLD NAPGEIRLKV DNLCGPGVND VSFVLRKGEI LGISGLMGAG RTELKMKVLYG AMPRTSGYVT LDGHEVVTRS PQDGLANGIV YISEDKRKRDG LVLGMSVKEN MSLTALDYFS RAGGSLKHKD EQQAVGDFIR LFNVKTPSME QAIGLLSGGN QQKVAIARGL MTRPKVLILD EPTRGVDVGA KKEIYQLINQ FKADGLSIL VSSEMPEVLG MSDRIIVMHE GHLSGEFTRE QATQEVLMMAA AVGKLN RVNQ E
<b>Source</b>	E.coli
<b>Target Names</b>	rbsA
<b>Protein Names</b>	Recommended name: Ribose import ATP-binding protein RbsA EC= 3.6.3.17
<b>Expression Region</b>	1-501
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.