



# Recombinant Salmonella paratyphi A UDP-N-acetylglucosamine 1-carboxyvinyltransferase (murA)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP722137SBM                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Abbreviation</b>      | murA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                   |
| <b>Uniprot No.</b>       | Q5PLD1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Immunogen Species</b> | Salmonella paratyphi A (strain ATCC 9150 / SARB42)                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Sequence</b>          | MDKFRVQGPT TLQGEVTISG AKNAALPILF AALLAEEPVE IQNVPKLDV<br>DTSMKLLSQL GAKVERNCSV HIDASQVNVF CAPYDLVKTM RASIWALGPL<br>VARFGQGQVS LPGGCTIGAR PVDLHITGLE QLGATIKLEE GYVKASVEGR<br>LKGAHIVMDK VSVGATVTIM CAATLAEGTT IIENAAREPE IVDTANFLVT<br>LGAKIAGQGT DRITIEGVER LGGGVYRVLP DRIETGTFLV AAAISRGKIL<br>CRNAQPDTLD AVLAKLRDAG ADIEVGEDWI SLDMHGKRPK AVNVRTAPHP<br>AFPTDMAQAF TLLNLVAEGT GFITETVFEN RFMHVPELSR MGARAEIESN<br>TVICHGVETL SGAQVMATDL RASASLVLAG CIAEGTTIVD RIYHIDRGYE<br>RIEDKLRLALG ANIERVKGE |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Target Names</b>      | murA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Protein Names</b>     | Recommended name: UDP-N-acetylglucosamine 1-carboxyvinyltransferase<br>EC= 2.5.1.7 Alternative name(s): Enoylpyruvate transferase UDP-N-<br>acetylglucosamine enolpyruvyl transferase Short name= EPT                                                                                                                                                                                                                                                                                                |
| <b>Expression Region</b> | 1-419                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                          |



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.