



Recombinant Idiogrammarina loihensis 50S ribosomal protein L4 (rpLD)

Product Code	CSB-MP722219IAAA
Abbreviation	rpLD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5QXZ1
Product Type	Recombinant Protein
Immunogen Species	Idiogrammarina loihensis (strain ATCC BAA-735 / DSM 15497 / L2-TR)
Purity	≥85% (SDS-PAGE)
Sequence	MELTLKDAKG ALEVSEATFG REFNEALVHQ VVVAYAAGAR QGTKAQKTRS EVAGGGKKPW RQKGTGRARA GTIRSPIWRS GGATFAAKPQ NHSQKVNKKM YRGAIKSILS ELIRQERLIV VEKFGVDEPK TKQLAAKLKE MDLNDVLIVT KEVDENLFLA SRNLHKVDVR DVQGIDPVSL IAFEKVLMTA DAVKQLEEV L S
Source	Mammalian cell
Target Names	rpLD
Protein Names	Recommended name: 50S ribosomal protein L4
Expression Region	1-201
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.