



Recombinant Pongo abelii Tetratricopeptide repeat protein 35 (TTC35)

Product Code	CSB-EP722333PYX-B
Abbreviation	TTC35
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R882
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	AKVSELYDV TWEEMRDKMR KWREENSRNS EQIVEVGEEL INEYGSKLGD DIWIIYEQVM IAALDYGRDD LALFCLQELR RQFPGSHRVK RLTGMRFEAM ERYDDAIQLY DRILQEDPTN TAARKRKIAI RKAQGKNVEA IRELNEYLEQ FVG DQEA WHE LAELYINEHD YAKAAFCLEE LMMTNPHNHL YCQQYAEVKY TQGGLENLEL SRKYFAQALK LNNRNMRALF GLYMSASHIA SNPKASAKTK KDNMKYASWA ASQINRAYQF AGRSKKETKY SLKAVEDMLE TLQITQS
Source	E.coli
Target Names	EMC2
Protein Names	Recommended name: Tetratricopeptide repeat protein 35 Short name= TPR repeat protein 35
Expression Region	2-297
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.