



# Recombinant Pongo abelii Endoplasmic reticulum lectin 1 (ERLEC1)

<b>Product Code</b>	CSB-MP722344PYX
<b>Abbreviation</b>	ERLEC1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5R8S4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	LPQLSDD IPFRVNWPGT EFSLPTTGVL YKEDNYVIMT TAHKEKYKCI LPLVTSGDEE EEKDYKGNP RELLEPLFKQ SCSYRIESY WTYEVCHGKH IRQYHEEKET GQKINIHEY Y LGNMLAKNLL FEKEREAEK EKSNEIPTKN IEQMTPYYP VGMGNGTPCS LKQNRPRSST VMYICHPEK HEILSVAEVT TCEYEVVILT PLLCSHPKYR FRASPVNDIF CQSLPGSPFK PLTLRQLEQQ EEILRVPFRR NKEEDLQSTK EERFPaihkp IAIGSQPVLT VGTTHISKLT DDQLIKEFLS GSYCFHGGVG WWKYEFcygk HVHqYHEDKD SGKTSVVVGT WNQEEHIEWA KKNTARAYHL QDDGTQTVRM VSHFYGNNGDI CDITDKPRQV TVKLKCKESD SPHAVTVYML EPHSCQYILG VESPVICKIL DTADENGLLS LPN
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ERLEC1
<b>Protein Names</b>	Recommended name: Endoplasmic reticulum lectin 1 Alternative name(s): ER lectin Short name= Erlectin
<b>Expression Region</b>	34-483
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.