



Recombinant Pongo abelii Eukaryotic translation initiation factor 1A, X-chromosomal (EIF1AX)

Product Code	CSB-BP722369PYX
Abbreviation	EIF1AX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RA42
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	MPKNKGKGGK NRRRGKNENE SEKRELVFKE DGQEYAQVIK MLGNRLEAM CFDGVKRLCH IRGKLRRKVV INTSDIILVG LRDYQDNKAD VILKYNAD EA RSLKAYGELP EHAKINETDT FGPGDDDEIQ FDDIGDDDED IDDI
Source	Baculovirus
Target Names	EIF1AX
Protein Names	Recommended name: Eukaryotic translation initiation factor 1A, X-chromosomal Short name= eIF-1A X isoform Alternative name(s): Eukaryotic translation initiation factor 4C Short name= eIF-4C
Expression Region	1-144
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.