



# Recombinant Pongo abelii Pre-mRNA-splicing factor RBM22 (RBM22)

<b>Product Code</b>	CSB-MP722388PYX
<b>Abbreviation</b>	RBM22
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5RAY5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	ATSLGSNTY NRQNWEDADF PILCQTCLGE NPYIRVTKEK YGKECKICAR PFTVFRWCPG VRMRFKKTEV CQTC SKLKNV CQTCLLDLEY GLPIQVRDAG LSFKDDMPKS DVNKEYYTQN MEREISNSDG TRPVGMLGKA TSTSDMLLKL ARTTPYYKRN RPHICSFVVK GECKRGECEP YRHEKPTDPD DPLADQNIKD RYYGINDPVA DKLLKRASTM PRLDPPEDKT ITTLYVGGLG DTITETDLRN HFYQFGEIRT ITVVQRQQCA FIQFATRQAA EVAAEKSFNK LIVNGRRLNV KWGRSQAARG KEKEKDGTTD SGIKLEPVPVPG LPGALPPPPA AEEEASANYF NLPPSGPPAV VNIALPPPPG IAPPPPPGFG PHMFHPMGPP PPFMRAPGPI HYPSQDPQRM GAHAGKHSSP
<b>Source</b>	Mammalian cell
<b>Target Names</b>	RBM22
<b>Protein Names</b>	Recommended name: Pre-mRNA-splicing factor RBM22 Alternative name(s): RNA-binding motif protein 22
<b>Expression Region</b>	2-420
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.