



# Recombinant Pongo abelii Sestrin-2 (SESN2)

<b>Product Code</b>	CSB-YP722416PYX
<b>Abbreviation</b>	SESN2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5RCB4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MIVADSECRA ELKDYLRFP GGVGDSGPGE EQRESRARRG PRGPSAFIPV EEVLREGAES LEQHLGLEAL MSSGRVDNLA VVMGLHPDYF TSWRLHYLP LHTDGPLASS WRHYAIMAA ARHQCSYLVG SHMAEFLQTG GDPEWLLGLH RAPEKLRKLS EINKLLAHRP WLITKEHIQA LLKTGEHTWS LAELIQALVL LTHCHSLSSF VFGCGILPEG DADGSPAPQA PTPPSEQSSP PSRDPLNNSG GFESARDVEA LMERMRLQE SLLRDEGTSQ EEMESRFELE KSESLLVTPS ADILEPSHP DMLCFVEDPA FGYEDFTRRG AQAPPTFRAQ DYTWEDHGYS LIQRLYPEGG QLLDEKFQAA YSLTYNTIAM LSGVDTSVLR RAIWNYIHCV FGIRYDDYDY GEVNQLLERN LKVYIKTVAC YPEKTRRMV NLFWRHFRHS EKVHVNLLLL EARMQAALLY ALRAITRYMT
<b>Source</b>	Yeast
<b>Target Names</b>	SESN2
<b>Protein Names</b>	Recommended name: Sestrin-2
<b>Expression Region</b>	1-480
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.