



Recombinant Danio rerio WASH complex subunit CCDC53 (ccdc53)

Product Code	CSB-EP722479DIL-B
Abbreviation	ccdc53
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RGJ6
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MDADGLPIVG SGVDLTKVPA IQRRRIVAFL NQFIVHTVRF LNRFSSTVCEE KLSTVSLRIQ QIETTLSILE AKLASIPGLE EVTVDGVRAP AETNGPAADT SRATAPPAEA SQQTQAAPQE PKTEAPAENI MTVAKDPRYA RYLKMOVQGV PVMAIKNKMM MEGLDPNLLD TPDAAVPDAS KKRLEEQDDD SSGSESSFSD
Source	E.coli
Target Names	washc3
Protein Names	Recommended name: WASH complex subunit CCDC53 Alternative name(s): Coiled-coil domain-containing protein 53
Expression Region	1-200
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.