



Recombinant *Thermus thermophilus* GTPase obg (obg)

Product Code	CSB-YP722512TNT
Abbreviation	obg
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5SHE9
Product Type	Recombinant Protein
Immunogen Species	<i>Thermus thermophilus</i> (strain HB8 / ATCC 27634 / DSM 579)
Purity	≥85% (SDS-PAGE)
Sequence	MFQDVLVITV AAGRGGDGAV SFRREKFVPK GGPDGGDGGR GGSVYLRARG SVDSLSRLSK RTYKAEDGEH GRGSQQHGGR GEDLVIEVPR GTRVFDADTG ELLADLTEEG QTVLVARGGA GGRGNMHFVS PTRQAPRFAE AGEEGEKRRRL RLELMLIADV GLVGYPNAGK SSLLAAMTRA HPKIAPYPFT TLSPNLGVVE VSEERFTLA DIPGIEGAS EGKGLGLEFL RHIARTRVLL YVLDAADEPL KTLETLRKEV GAYDPALLRR PSLVALNKVD LLEEEAVKAL ADALAREGLA VLPVSALTGA GLPALKEALH ALVRSTPPPE MPKPVPKRKEV QAGVEVVPVA EGVYEVRAPE VERYLARIKG DLMEAAGYLQ EVFRRQGVEA ALRAKGVRAG DLVRIGGLEF EYIPEV
Source	Yeast
Target Names	obg
Protein Names	Recommended name: GTPase obg Alternative name(s): GTP-binding protein obg TT1381
Expression Region	1-416
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.