



Recombinant Mouse Histone deacetylase complex subunit SAP30L (Sap30l)

Product Code	CSB-MP722543MO
Abbreviation	Sap30l
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5SQF8
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MNGFSTEEDS REGPPAAPAA APGYGQSCCL IADGERCVRP AGNASFSKRV QKSISQKCLK LDIDKSVRHL YICDFHKNFI QSVRNKRKRK ASDDGGDSPE HDADIPEVDL FQLQVNTLRR YKRHYKLQTR PGFNKAQLAE TVSRHFRNIP VNEKETLAYF IYMVKSNRSR LDQKSEGSKQ LE
Source	Mammalian cell
Target Names	Sap30l
Protein Names	Recommended name: Histone deacetylase complex subunit SAP30L Alternative name(s): Sin3 corepressor complex subunit SAP30L Sin3-associated protein p30-like
Expression Region	1-182
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.