



Recombinant Mouse F-box/WD repeat-containing protein 11 (Fbxw11)

Product Code	CSB-EP722546MO
Abbreviation	Fbxw11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5SRY7
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MEPDSVIEDK TIELMCSVPR SLWLGCANLV ESMCALSCLQ SMPSVRCLQI SNGTSSVIVS RKRPESENYQ KEKDLCIKYF DQWSESDQVE FVEHLISRM HYQHGHINSY LKPMQLRDFI TALPEQGLDH IAENILSYLD ARSLCAAELV CKEWQRVISE GMLWKKLIER MVRTDPLWKG LSERRGWDQY LFKNRPTDGP PNSFYRSLYP KIIQDIETIE SNWRCGRHNL QRIQCRSENS KGVYCLQYDD DKIIISGLRDN SIKIWDKSSL ECLKVLTGHT GSVLCLQYDE RVIVTGSSDS TVRVWDVNTG EVLNTLIHNN EAVLHLRFSN GLMVTCSKDR SIAVWDMASA TDITLRRVLV GHRAAVNVVD FDDKYIVSAS GDRTIKVWST STCEFVRTLN GHKRGIAQLQ YRDRLVVSGS SDNTIRLWDI ECGACLRVLE GHEELVRCIR FDNKRIVSGA YDGKIKVWDL QAALDPRAPA STLCLRTLVE HSGRVFRLQF DEFQIISSSH DDTILIWDFL NVPPSAQNET RSPSRITYTYI SR
Source	E.coli
Target Names	Fbxw11
Protein Names	Recommended name: F-box/WD repeat-containing protein 11 Alternative name(s): F-box and WD repeats protein beta-TrCP2 F-box/WD repeat-containing protein 1B Homologous to Slimb protein Short name= HOS
Expression Region	1-542
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a member of the F-box protein family which is characterized by an approximately 40 amino acid motif, the F-box. The F-box proteins constitute one of the four subunits of ubiquitin protein ligase complex called SCFs (SKP1-cullin-F-box), which function in phosphorylation-dependent ubiquitination. The F-box proteins are divided into 3 classes: Fbws containing WD-40 domains, Fbls containing leucine-rich repeats, and Fbxs containing



either different protein-protein interaction modules or no recognizable motifs. This protein belongs to the Fbws class and, in addition to an F-box, contains multiple WD40 repeats. This gene contains at least 14 exons, and its alternative splicing generates 3 transcript variants diverging at the presence/absence of two alternate exons.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.