



# Recombinant Human Integrator complex subunit 11 (INTS11)

<b>Product Code</b>	CSB-MP722574HU
<b>Abbreviation</b>	INTS11
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5TA45
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MPEIRVTPLG AGQDVGRSCI LVSIAGKNVM LDCGMHMGFN DDRRFPDFSY ITQNGRLTDF LDCVIISHFH LDHCGALPYF SEMVGYDGPI YMTHPTQAIC PILLEDYRKI AVDKKGEANF FTSQMIKDCM KKVAVVHLHQ TVQVDDELEI KAYYAGHVLG AAMFQIKVGS ESVVYTG DYN MTPDRHLGAA WIDKCRPNLL ITESTYATTI RDSKRCRERD FLKKVHETVE RGGKVLIPVF ALGRAQELCI LLETFWERMN LKVPIYFSTG LTEKANHYK LFIPWTNQKI RKTFVQRNMF EFKHIKAFDR AFADNPGPMV VFATPGMLHA GQSLQIFRKW AGNEKNMVIM PGYCVQGTVG HKILSGQRKL EMEGRQVLEV KMQVEYMSFS AHADAKGIMQ LVGQAEPEPV LLVHGAEAKM EFLKQKIEQE LRVNCYMPAN GETVTLPTSP SIPVGISLGL LKREMAQGLL PEAKKPRLLH GTLIMKDSNF RLVSSEQALK ELGLAEHQLR FTCRVHLHDT RKEQETALRV YSHLKSVLKD HCVQHLPDGS VTVESVLLQA AAPSEDPGK VLLVSWTYQD EELGSFLTSL LKKGLPQAPS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	INTS11
<b>Protein Names</b>	Recommended name: Integrator complex subunit 11 Short name= Int11 EC= 3.1.27.-Alternative name(s): Cleavage and polyadenylation-specific factor 3-like protein Short name= CPSF3-like protein Protein related to CPSF subunits of 68
<b>Expression Region</b>	1-600
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.