



Recombinant Human Zinc finger and BTB domain-containing protein 37 (ZBTB37)

Product Code	CSB-EP722578HU-B
Abbreviation	ZBTB37
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5TC79
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEKGGNIQLE IPDFSNSVLS HLNQLRMQGR LCDIVNVVQG QAFRAHKVVV AASSPYFRDH MSLNEMSTVS ISVIKNPTVF EQLLSFCYTG RICLQLADII SYLTAASFLQ MQHIIDKCTQ ILEGIHFKIN VAEVEAELSQ TRTKHQERPP ESHRVTPNLN RLSLSPRHNTK KGNRRGQVSA VLDIRELSPP EESTSPQIIE PSSDVESREP ILRINRAGQW YVETGVADRG GRSDDEVRLV GAVHIKTENL EEWLGPENQP SGEDGSSAEE VTAMVIDTTG HGSVGQENYT LGSSGAKVAR PTSSEVDRFS PSGSVVPLTE RHRARSESPG RMDEPKQPSS QVEESAMMGV SGYVEYLREQ EVSERWFRYN PRLTCIYCAK SFNQKGSLLDR HMRLHMGITP FVCRMCGKKY TRKDQLEYHI RKHTGNKPFH CHVCGKSFPF QAILNQHFRK NHPGCIPLEG PHSISPETTV TSRGQAEES PSQEETVAPG EAVQGSVSTT GPD
Source	E.coli
Target Names	ZBTB37
Protein Names	Recommended name: Zinc finger and BTB domain-containing protein 37
Expression Region	1-503
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.