



Recombinant Danio rerio Sel1 repeat-containing protein 1 (selrc1)

Product Code	CSB-EP722592DIL-B
Abbreviation	selrc1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5TYQ3
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MAGLINFED E QEVKQFLDNL GVEYSYQCYR EKDPEGCHRL ADYLEGVKKN YESTAQVLQH NCEVNAHAQS CYKLGAYHVT GKGGMKKCLK TAYSCFLKSC NTQGKKSVD A CHNVGLLAQD GRALETGPD T VARQYFEKA CEGGFAPSCF NLSTLYIQGF PGLDKSMPLA LKYALKACDL GHVWGCANAS RMYKLG DGTD KDEQRAEELK NRAKDLHGQE KERQLKFGE
Source	E.coli
Target Names	coa7
Protein Names	Recommended name: Sel1 repeat-containing protein 1 Alternative name(s): Beta-lactamase hcp-like protein
Expression Region	1-229
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.