



Recombinant Rat Tetratricopeptide repeat protein 26 (Ttc26)

Product Code	CSB-BP722599RA
Abbreviation	Ttc26
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5U2N8
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	MMLSRAKPAV GGESPHTDKR KKKGRKIPKL EDLLSQRDFT GAITLLEFKR HVGEQEEDTN LWIGYCAFHL GDYKRALEEY ENAAKEENCN PEVWVNLACT YFFLGMYKQA EAAGFKAPKS RLQNRLLFHL AHKFNDEKKL MNFHQNLQDI KEDQLSLASI HYMRSHYQEA IDIYKRILLD NREYLALNVY VALCYKLDY YDVSQEV LAV YLQQIPDSTI ALNLKACNHF RLYNGKAAEA ELKSLMDNAS SPFEFAKELI RHNLVFRGG EGALQVLPPL VDVPEARLN LVIYYLRQDD VQEAYNLIKD LEPTTPQEYI LKGVVNAALG QEMGSRDHMK IAQQFFQLVG GSASECDTIP GRQCMASCFE LKQFDDVLI YLNSFKSYFY NDDIFNFNYA QAKAATGNTS EGEEVLLIQ SEKLKNDYIY LSWLARCIM NKKPRLAWEL YLMETSGES FSLLQLIAND CYKMGQFYYS AKAFDVLRL DPNPEYWEGK RGACVGFQIM ILAGREPKET LREVLHLLRS TGNTQVEYII RIMKKWAKEN RVPI
Source	Baculovirus
Target Names	Ttc26
Protein Names	Recommended name: Tetratricopeptide repeat protein 26 Short name= TPR repeat protein 26
Expression Region	1-554
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.