



# Recombinant Rat Tetratricopeptide repeat protein 26 (Ttc26)

<b>Product Code</b>	CSB-MP722599RA
<b>Abbreviation</b>	Ttc26
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5U2N8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MMLSRAKPAV GGESPHTDKR KKKGRKIPKL EDLLSQRDFT GAITLLEFKR HVGEQEEDTN LWIGYCAFHL GDYKRALEEY ENAAKEENCN PEVWVNLACT YFFLGMYKQA EAAGFKAPKS RLQNRLLFHL AHKFNDEKKL MNFHQNLQDI KEDQLSLASI HYMRSHYQEA IDIYKRILLD NREYLALNVY VALCYKLDY YDVSQEV LAV YLQQIPDSTI ALNLKACNHF RLYNGKAAEA ELKSLMDNAS SPFEFAKELI RHNLVFRGG EGALQVLPPL VDVPEARLN LVIYYLRQDD VQEAYNLIKD LEPTTPQEYI LKGVVNAALG QEMGSRDHMK IAQQFFQLVG GSASECDTIP GRQCMASCFE LKQFDDVLI YLNSFKSYFY NDDIFNFNYA QAKAATGNTS EGEEVLLIQ SEKLKNDYIY LSWLARCIM NKKPRLAWEL YKMETSGES FSLLQLIAND CYKMGQFYYS AKAFDVLRL DPNPEYWEGK RGACVGFQIM ILAGREPKET LREVLHLLRS TGNTQVEYII RIMKKWAKEN RVPI
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Ttc26
<b>Protein Names</b>	Recommended name: Tetratricopeptide repeat protein 26 Short name= TPR repeat protein 26
<b>Expression Region</b>	1-554
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.