



Recombinant Rat Lebercilin (Lca5)

Product Code	CSB-MP722606RA
Abbreviation	Lca5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5U2Y9
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	<p> MGERDQSSDT DQERKGGKHP YSSHYSSDFG SSPQSSGPSS PVNTPPCASN REKNPKRHLS DHHVRHQTVR KVSPKAIPNK KGVRVGFERSQ SLNREPLRKD PDIVTKRVLS ARLLKITELQ NEVSELQVKL AQLLKENKAL KSLQYRQEKA LSKFEDAENE ISQLIHRHNG EITALKERLR KSQEKERMAE KRVKDTEGEL FRTKFSLQKL KKISEARHLP ERDDLAKKLV SAELKLDNTE RKIKELSKNL ELSTNSFQRQ LLAERKRAFE AHDENKILQK ELQRLHHKLIK EKEKELDIKN IYANRLPKSS PKKEKEIARK HASCQSDFTD LCTKAVQTAE DFEPEEFPFT AQTVLCYENR WDGSEYLPSY LEYQERNEHG SESLSSTLEQ EGKYSDDQDS CTAKQETGKS ECEWEREELD KVKGKSSLLG RTEKLVLEAG SVQTESYQAQ TIDKFQDEAE RLKTEMLLAK LNEINKELQD PQSLSRPLLP LLPNFESKLY SPDRSLRPYT FSESLDRSFN GQHLSDMSFL TPRDEGQSPG LIRSPAPSDE FSGSYVPSF GKTPGKSNPP SQKNSLLDFQ SNMSEGLSKD SLDFISRKEK KATLMEQLFG TSASASATNT SASTTTNTTT NTNTNTNTNT NTNCNVPSKS IDSHFPAASR GDLPLHFLS GDRSSRVREP DGEDEDLFLS EGRSFNPNRH RLKHTSTKPT VTAVDSIDED IEEVTLR </p>
Source	Mammalian cell
Target Names	Lca5
Protein Names	Recommended name: Lebercilin Alternative name(s): Leber congenital amaurosis 5 protein homolog
Expression Region	1-727
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a protein that is thought to be involved in centrosomal or ciliary functions. Mutations in this gene cause Leber congenital amaurosis type V. Alternatively spliced transcript variants are described.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.