



Recombinant Human Uncharacterized protein C12orf34 (C12orf34)

Product Code	CSB-MP722624HU
Abbreviation	C12orf34
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5U5X8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MLACLQRTQN APGQHLACPS KSLELRKCEA VASAMHSSRY PSPAELDAYA EKVANSPLSI KIFPTNIRVP QHKHLSRTVN GYDTSGQRYS PYPQHTAGYQ GLLAIVKAAV SSSSTAAPAG PAKSVLKSAE GKRTKLSPAA VQVGIAPYPV PSTLGPLAYP KPPEAPAPPP GLPAAATAAS VIPLPGRGLP LPPSNLPSIH SLLYQLNQQC QAPGAAPPAC QGMAIPHPSP AKHGVPVSFP SMAYSAAAGL PDCRKGTELG QGATQALTLA GAAKPAGYAD SGLDYLLWPQ KPPPPPPQPL RAYSGSTVAS KSPEACGGRA YERASGSPLN CGVGLPTSFT VGQYFAAPWN SVLVTPTSDC YNPAAAVVVT ELGPGAAREL AGPPADALSG LPSKSVCNTS VLSSSLQSLE YLINDIRPPC IKEQMLGKGY ETVAVPRLLD HQHAHIRLPV YR
Source	Mammalian cell
Target Names	FAM222A
Protein Names	Recommended name: Uncharacterized protein C12orf34
Expression Region	1-452
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.