



Recombinant Acanthamoeba polyphaga mimivirus

Uncharacterized protein R155 (MIMI_R155)

Product Code	CSB-MP722669ADAZ
Abbreviation	MIMI_R155
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UPL5
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MYSFYKPFGE VNTDINKPID KPIMSCSTHT YNNDKLVYKQ FKRLTKLNND LKKFCDDPDF VNSTFRRNII YNIHNKRIHL IKKLSMIEPD ITYNKFIIMW TKNSFNKIKS TSRKVPKEAK SDWYLYKNFV FVQGLNRHKQ TV
Source	Mammalian cell
Target Names	MIMI_R155
Protein Names	Recommended name: Uncharacterized protein R155
Expression Region	1-142
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.