



Recombinant Acanthamoeba polyphaga mimivirus

Uncharacterized protein R325 (MIMI_R325)

Product Code	CSB-BP722717ADAZ
Abbreviation	MIMI_R325
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UQR6
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	SYVKSDITNK SYLVQNLPNK EEEAHILGII DKRISILDDY LRENINKFTE YKPYIEQFNN RIKKLTLYEN APDGKYTSFT VDKGKEIALC LRSRKTGQIH DINLVMYVTL HELAHVACPE TDHTELFKKI FIFLIKISID LNIYKHIDYE ADPAEYCGLV IDEDLLK
Source	Baculovirus
Target Names	MIMI_R325
Protein Names	Recommended name: Uncharacterized protein R325
Expression Region	31-197
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.