



Recombinant Acanthamoeba polyphaga mimivirus Uncharacterized protein R650 (MIMI_R650)

Product Code	CSB-BP722727ADAZ
Abbreviation	MIMI_R650
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UR07
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MIQNNPKSIG SSSNKSARSS GSRQNILKMI QFENQRLHEY YNDFSSLKET ILTAPSIVEY LFQKPSMKYD TIEKIYKLFT DPKSNYLENY TNITRVTKYI SDYLFADMIE LIELQDKYHV CNDNYNYVVD MILKNRNKYP KYTYLSQRYT PLIQTIGDTL DNIQLVEKYV DDDGAKISIK YEISWIYNWT FKNDCKYNLD IPFVDFDFVV LVYRERLILF VIEVDRDINN PLYHTESIHK QILFQMNVF VLRISQQHLK DNLKKRIIRF IKNVSSKVY LSEKYITPNI NLIDMDLVKL ELEKFYQSYK CCNSYYHKSL LDNSTIYPDP IKNKFHNGRN MGITLKTDKI YNKDDDEYFD TLIGTKFTRE PDPIVVTDD VFDWIINKCP KEPSKKQYSE RKISEKQQKY NDICISLVGK II
Source	Baculovirus
Target Names	MIMI_R650
Protein Names	Recommended name: Uncharacterized protein R650
Expression Region	1-422
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.