



Recombinant Haloarcula marismortui 1- (5-phosphoribosyl)-5-[(5-phosphoribosylamino)methylideneamino] imidazole-4-carboxamide isomerase (hisA)

Product Code	CSB-MP722754HTJ
Abbreviation	hisA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UZH1
Product Type	Recombinant Protein
Immunogen Species	Haloarcula marismortui (strain ATCC 43049 / DSM 3752 / JCM 8966 / VKM B-1809) (Halobacterium marismortui)
Purity	≥85% (SDS-PAGE)
Sequence	MYPEFEVPPA VDMQDGQVVQ LVGGERGTEK TYGDPVEAAQ RWVDAGARTL HLVDLDGAFE GERQNAAAID AVLDAVGADV DVQLGGGIRT AEDAVSLLDR GLDRVILGTA AVETPEIVGE ISDEHPGSVL VSLDAKDGEV VVSGWTEGTG LDPADAAGRY ADLGAGGILF TDVDVEGQLD GVRTDPVRRRL VDSVDIPVIA SGGVATINDV LALRSAGAAA VVVGSALEYEG QFTLDAAIDA LGESSE
Source	Mammalian cell
Target Names	hisA
Protein Names	Recommended name: 1-(5-phosphoribosyl)-5-[(5-phosphoribosylamino)methylideneamino] imidazole-4-carboxamide isomerase EC= 5.3.1.16 Alternative name(s): Phosphoribosylformimino-5-aminoimidazole carboxamide ribotide isomerase
Expression Region	1-246
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.